



[www.wheelchair2u.my](http://www.wheelchair2u.my)



# **PRODUCT CATALOGUE**



Founded in 2015, Good Take Care Trade registered with Kementerian Kewangan Malaysia (MOF) & Suruhanjaya Syarikat Malaysia (SSM) based in Malaysia.

The origin word of Good & Take Care primarily means

“**Satisfactory  
& Look after**”

Good & Take Care is our philosophy that emphasizes all of our products.

As we know, not good mobility can be challenged for the disabled person in daily life but good products of mobility can be easier and helpful to go through the whole day.

## ***Founder story.***

Everything started in 2015, the founder of Gospeedi losing mobility after diagnosed with Spinal Cord Injury at 42 years old. Becoming wheelchair-bound is a scary life challenge but with support from family and the spreading of positive energy, the founder successfully adopted a new way of life and overcome challenges.

During the period of adopting a new way of life, the founder observing something in his mind. He chooses to make custom wheelchair with a combination of the selected component such as different front custer, rear rim, sportrim and type of tire to make more practical and stylish.

As an active mountain bike bicycle before paralyzed, the founder looking for an idea of how group of strong hand paraplegia and another group of elderly people to move around faster, classy, well balanced with normal society.

And with all his knowledge, he put an effort to bring satisfactory products and fitting clients needs to lead them in fulfilling life regardless of the disability.

## Our Vision

Dedicated to continual improvement and contribute good disability products in Malaysia.

## Our Mission

Improve the quality life of disable people and senior citizen in Malaysia.

## Our Value

We believe combination of hard work and positive attitude will bring smile to our client.

## Trusted Clients.



# Gospeedi Motor.



Attach to chair in just second without transfer.



## 2 Options

- Front tire (Size 12 inch)
- Front tire (Size 16 inch)



## Lithium Battery

- Capacity 36v 8ah
- Airline requirement approved





### Smart Controller Display

- Button ON and OFF
- 3 speed & 1 reverse gear
- Display gear , time , total KM range , battery gauge



### Braking System

- Double disc brake
- Push button for handbrake
- Left & right handle system safety



### Body Frame

- Snap on & Snap off
- Adjustable height
- Safety connector to wheelchair
- Adjustable claw for different wheelchair tube



**“Freedom to move  
& Regain your independence”**

# Spareparts & Accesories

## Gospeedi Motor



Break Clip



Connector



Half Trottle



Smart Controller  
Display



Main Wire Harness



Battery 36v 8ah



Body Clam



Tire 12 inch



Tire 16 inch



Mud Guard



Stand Wheel



LED Front Lamp



Horn & Lamp Cable



Safety Switch



Charger 36v

# Gospeedi Rigid Frame Wheelchair.

**CustomMade**



## Standard

Standard wheel tire 24 inch



## + Gospeedi Wheelz

Wheelz 12 spokes with 7 color  
Schwalbe Right Run tire 24inch



## + Gospeedi Wheelz + Offroad Tire

Wheelz 12 spokes with Black color  
Branded tire 24 x 1.95 inch



## + Fat Bob Tire

Schwalbe fat tire 24 x 3 inch

## Product Details

- Available width size **14 and 16 inch**
- Front castor 4 inch
- Detachable side guard
- Quick-release rear wheel
- Material Aluminium Alloy
- Scissors lock
- Armrest
- 2 inch Cushion
- Body Cover
- Foldable Backrest



# Gospeedi Folding Frame Wheelchair.

CustomMade



## Standard

Standard wheel tire 24 inch



## + Gospeedi Wheelz

Wheelz 12 spokes with 7 color  
Schwalbe Right Run tire 24inch



## + Gospeedi Wheelz + Offroad Tire

Wheelz 12 spokes with Black color  
Branded tire 24 x 1.95 inch



## + Fat Bob Tire

Schwalbe fat tire 24 x 3 inch

## Product Details

- Available width size **14, 16 & 18 inch**
- Front castor 4 inch with Red fork
- Quick-release rear wheel
- Material Aluminium Alloy
- Push to lock
- Body Cover
- 2 inch Cushion
- Armrest

# Spareparts & Accesories

## Manual Wheelchair



Tire Schwable  
RightRun 24 inch



Tire IVStone 24 Inch  
with Tube



Front Caster  
4 inch & 5 inch




Sub Quick Release  
105, 115 & 120 mm




Sub Quick Release  
Sleeve



Bearing Rear Wheel  
6201z



Bearing Caster  
608 2rs



Bearing Fork  
6001 2rs



Bearing Rear Wheel  
R8 2rs



Fork Double Bearing



Fork



Fork with Stem



Fork Sub Screw



Rigid Footrest



Folding Footrest



Scissors Brake  
Type 1



Scissors Brake  
Type 2



Push Brake Type 1



Push Brake Type 2



Wheel Spoke Cover



Cover Body  
Shor 10.5 inch



Cover Body  
Long 15 inch



Cover Body  
Short 11 inch



Basket Nett



Stokeen



Foot Guard  
2 x 24 inch



Foot Guard  
4 x 30 inch

# Gospeedi Electric Wheelchair.

## Remote Control



### Intelligent Controller

- 5 speed gear
- Fold & Unfold button
- 360 degree controller joystick

### Luggage Style

- Easy storage & carry
- Fit into most car trunk



### 500w Motor Power

- Durable motor
- Switchable Modes  
Auto to Manual




### Remote Controller

- Control wheelchair wireless
- Auto folding & unfolding



### Battery Capacity

- 8 hours full charging
- 24v 12ah Lithium
- Display gear , time , total KM range , battery gauge
- Airline requirement approved



### Shock Absorber

- Comfortable ride
- Push button for handbrake
- Left & right handle system safety



### Carry Organizer

- Side Storage Bag
- Tools
- Charger



“ *Comfortable ride venture off pavement & outdoor asphalt* ”

# Spareparts & Accesories

## Electric Wheelchair



Intelligent Controller



250w Motor



LED Light with  
USB Port



Front Tire with Bearing  
(8 x 2 inch)



Rear Tire with Tube  
(12.5 x 1/4 inch)



Battery 24v 12ah

# Gospeedi Wheelz.



## Gospeedi Wheelz

- Size 24 inch
- 12 spokes
- 7 stylish color
- Schwalbe Right Run tire 24 inch with tube



## + Gospeedi Offroad Tire

- Size 24 inch
- Black spokes color
- Aluminum anodized handrim
- Kenda tire 24 x 1.95 inch with tube



## + Gospeedi Fat Tire

- Size 24 x 3 inch
- Schwalbe fat tire 24 inch with tube
- Black aluminum anodized handrim
- Quick Release button and Spacer

# Multifunctional Chair.



360 degree  
Backrest Opening




Double Lock  
• Safety belt  
• Fix setting position



Height Adjustment  
• 9 steps adjustment  
• Manual lift








### Silent Wheel

- Easy control
- Foot braking for safety



### Waterproof Material

- Medical grade
- Molded sitting
- Soft cushion backrest

**“It’s a Wheelchair, Transfer chair, Commode chair and Shower”**

# ROHO Cushion.

The High and Low-Profile Cushion by ROHO effectively redistribute pressure and conform to the user's body in order to prevent ulcers and tissue breakdown.



It's soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated.



## ***Please choose your suitable level***

- Available size 14 and 16 inch
- Low and High Profile

# Gospeedi Exclusive T-Shirt



Size Available | XS | S | M | L | XL | XXL

## ***T-Shirt Details***

- Microfiber material
- Activewear
- Quick dry
- CoolFit



**GOOD TAKECARE TRADE** (NSD192746H)

No . 37, Jalan BM 4/3, Seksyen 4, Bandar Bukit Mahkota, 43000 Kajang, Selangor



**+6019 224 7300**



**motorwheelchair**



**[www.wheelchair2u.my](http://www.wheelchair2u.my)**